


# MoxieDori Face Mask Extender



## INSTRUCTIONS

Turn your extender so the logo is right side up and it's "smiling"



Insert mask elastic from one side of the mask into the extender tension adjustment of your choice and gently tug to secure.



Put mask on, following standard safety protocol, with the extender still attached.



Secure remaining side mask elastic around the same extender tension adjustment choice.



Depending on the style and size of mask, adjust extender to where it comfortably rests on the head.



While your mask may be disposable, your extender can be cleaned, sterilized and is reusable!

Questions? Comments?  
[info@moxiedori.com](mailto:info@moxiedori.com)

**MOXIEDORI**